# NGAL Rapid ELISA Kit

EN

DE

FR

IT

ES

DA

SE

RU

**KIT 037**

**BIOPORTO**<sup>®</sup>  
Diagnostics

**Revision: NR2012-04EN****Please read these instructions carefully****INTENDED USE**

For the *in-vitro* determination of human NGAL in urine or plasma as a marker of acute kidney injury which may lead to acute renal failure.

**INTRODUCTION**

NGAL (neutrophil gelatinase-associated lipocalin) was first discovered in 1989<sup>1</sup> and, as its name implies, belongs to the lipocalin family of proteins. These are typically small secreted proteins characterized by their ability to bind hydrophobic molecules in a structurally conserved pocket formed by  $\beta$ -pleated sheet, to bind to specific cell-surface receptors, and to form macromolecular complexes. NGAL was fully characterized and named in 1993<sup>2</sup>, but has many synonyms: NL (neutrophil lipocalin; HNL: human NL)<sup>3</sup>, lipocalin 2, siderocalin, oncogene protein 24p3<sup>4</sup> or uterocalin<sup>5</sup> (in the mouse) and *neu*-related lipocalin<sup>6</sup> or 25-kDa  $\alpha_2$ -microglobulin-related protein<sup>7</sup> (in the rat). Human NGAL consists of a single disulfide-bridged polypeptide chain of 178 amino-acid residues with a calculated molecular mass of 21 kDa<sup>2</sup>, but glycosylation increases its apparent molecular mass to 25 kDa. In neutrophils (neutrophilic polymorphonuclear leukocytes) it occurs in monomer and homodimer, forms with a small percentage of higher molecular weight forms, and some of it is found in complex with 92-kDa human neutrophil type IV collagenase (gelatinase B or matrix metalloproteinase-9, MMP-9)<sup>8</sup>.

NGAL was originally isolated from the supernatant of activated human neutrophils<sup>1</sup>, but it is also expressed at a low level in other human tissues including the kidney, prostate and epithelia of the respiratory and alimentary tracts<sup>9,10</sup>. It is strongly expressed in adenomas and inflamed epithelia of the bowel<sup>11</sup>, adenocarcinomas of the breast<sup>12</sup>, and urothelial carcinomas<sup>13</sup>.

Because of its small molecular size and resistance to degradation, NGAL is readily excreted and detected in the urine, both in its free form and in complex with MMP-9. Urinary levels correlate with

plasma levels whatever the cause of increased NGAL production<sup>14</sup>, but particularly high urinary levels can be expected when NGAL is released directly into the urine by the kidney tubules or urothelial carcinomas. It is uncertain how far NGAL/MMP-9 complexes from sources remote from the urinary tract are excreted as such into the urine or reform in the urine after independent excretion of NGAL and MMP-9<sup>8</sup>.

While the functions of NGAL are not fully understood, NGAL appears to be up-regulated in cells under "stress", e.g. from infection, inflammation or neoplastic transformation, or in tissues undergoing involution, such as the postpartum mouse uterus and mammary glands on weaning<sup>5</sup>. In relation to a possible antibacterial role, it binds enterobactin and other siderophores, depriving the microorganisms of Fe<sup>3+</sup>, an important nutritional requirement<sup>15</sup>. Its complex formation with MMP-9 appears to protect MMP-9 enzymatic activity from degradation<sup>8</sup>. The up-regulation of NGAL in involuting tissues has led to the postulation of a role in apoptosis, but it appears more likely that NGAL is associated with a survival response<sup>16</sup>.

**NGAL and the kidney.** As early as 1989, NGAL (24p3) was found to be expressed by mouse kidney cells and to undergo an early, dramatic up-regulation (14- to 20-fold) in response to SV40 viral infection<sup>17</sup>. A similar early and dramatic up-regulation was later observed in other types of kidney injury, and raised plasma levels of NGAL were found to be strongly correlated with decreased renal function in patients with renal damage due to systemic vasculitis<sup>18</sup>. Urinary NGAL levels serve as an early marker of acute kidney injury after cardiopulmonary bypass surgery<sup>19,20</sup>, and both urinary and plasma levels of NGAL provide an early indication of acute renal injury in unselected patients in intensive care<sup>21</sup>. Raised urinary and plasma NGAL levels have also been observed in patients with chronic kidney diseases<sup>22</sup>, and patients with functioning renal grafts also showed raised urinary levels (detectable by Western blotting)<sup>13</sup>. It is therefore apparent that a large variety of renal disorders are associated with raised plasma and urinary levels of NGAL.

While plasma and urinary NGAL levels are closely correlated in acute conditions, it is to be expected that urinary NGAL levels will be particularly high after kidney injury severe enough to result in acute renal failure, acute tubular necrosis or acute tubulointerstitial nephropathy. However, the use of NGAL as a marker for renal disorders is subject to the proviso that other concurrent conditions that are independently associated with raised NGAL levels must be taken into account.

**NGAL in inflammation or infection.** NGAL is released from the specific granules of activated neutrophils<sup>1,2</sup> and serum levels rise in inflammatory or infective conditions, especially in bacterial infections<sup>23</sup>. Serum NGAL has thus been proposed as a marker of bacterial infection. However, as levels of NGAL may also be raised in neoplastic conditions and renal disorders independently of any infective process, this proposed application should be treated with caution. NGAL may also be raised in infections in patients with an uncountably low number of neutrophils due to leukemia or treated leukemia, showing that the source of the raised NGAL in infections is not only the neutrophils.

**NGAL and neoplasia.** The various types of cancer in which NGAL may be up-regulated (often with MMP-9) have been referred to above. This has been shown by its expression in tumor cells and its raised urinary levels, both in the free form and in complex with MMP-9<sup>8</sup>. Indeed, it has been proposed that urinary NGAL/MMP-9 complexes may serve as a marker of disease status for breast cancer patients<sup>24</sup>. Plasma levels have not usually been measured in these cases.

## PRINCIPLE OF THE ASSAY PROCEDURE

The assay is an ELISA performed in microwells coated with a monoclonal capture antibody to human NGAL. Bound NGAL is detected with a horseradish peroxidase (HRP)-conjugated monoclonal detection antibody and the assay is developed by incubation with a color-forming substrate. The assay is a rapid 2-step procedure:

**Step 1.** Aliquots of calibrators, diluted samples and any controls are incubated with HRP-conjugated

detection antibody in the coated microwells. Only NGAL will bind to both coat and detection antibody, while unbound materials are removed by washing.

**Step 2.** A chromogenic peroxidase substrate containing tetramethylbenzidine (TMB) is added to each test well. The HRP linked to the bound detection antibody reacts with the substrate to generate a colored product. The enzymatic reaction is stopped chemically, and the color intensity is read at 450 nm in an ELISA reader. The color intensity (absorbance) is a function of the concentration of NGAL originally added to each well. The results for the calibrators are used to construct a calibration curve from which the concentrations of NGAL in the test specimens are read.

## KIT COMPONENTS

Item	Contents	Quantity
①	12 x 8 coated Microwells + Frame	96 wells
②	5x Sample Diluent Conc.	1 x 60 mL
③-⑤	NGAL Rapid Calibrator 1-6 0, 0.2, 2, 5, 10, 20 ng/mL	6 x 1 mL
④	25x Wash Solution Conc.	1 x 30 mL
⑤	HRP-conjugated NGAL Antibody	1 x 6 mL
⑥	TMB Substrate	1 x 12 mL
⑦	Stop Solution	1 x 16 mL
⑧	Polypropylene U-Microwell Plate	96 wells

**Note:** Liquid reagents contain preservatives and may be harmful if ingested.

## MATERIALS REQUIRED BUT NOT PROVIDED

- Adjustable micropipettes covering the range 1-1000  $\mu$ L and corresponding disposable pipette tips
- Polypropylene tubes to contain up to 1000  $\mu$ L
- Tube racks
- Adjustable 8- or 12-channel micropipette (50-250  $\mu$ L range) or repeating micropipette (optional)

5. Clean 1 L and 500 mL graduated cylinders
  6. Deionized or distilled water
  7. Cover for microwell plate
  8. Clean container for diluted Wash Solution
  9. Apparatus for filling wells during washing procedure (optional)
  10. Lint-free paper towels or absorbent paper
  11. Disposable pipetting reservoirs
  12. Timer (60-minute range)
  13. Calibrated ELISA plate reader capable of reading at 450 nm (preferably subtracting reference values at 650 or 620 nm)
  14. Sodium hypochlorite (household bleach 1:10 dilution) for decontamination of specimens, reagents and materials
- contaminated specimens may give erroneous results.
10. Do not dilute clinical specimens directly in the coated microwells.
  11. Do not touch or scrape the bottom of the coated microwells when pipetting or aspirating fluid.
  12. Incubation times and temperatures other than those specified may give erroneous results.
  13. Do not allow the wells to dry once the assay has begun.
  14. The TMB Substrate is light sensitive. Keep away from bright light.
  15. Do not reuse microwells or pour reagents back into their bottles once dispensed.

## PRECAUTIONS

### For *in vitro* diagnostic use only

1. This kit should only be used by qualified laboratory staff.
2. Use separate pipette tips for each sample, calibrator and reagent to avoid cross-contamination.
3. Use separate reservoirs for each reagent. This applies especially to the TMB Substrate.
4. After use, decontaminate all specimens, reagents and materials by soaking for at least 30 minutes in sodium hypochlorite solution (household bleach diluted 1:10).
5. To avoid droplet formation during washing, aspirate the wash solution into a bottle containing bleach.
6. Avoid release into the environment. Dispose of containers and unused contents in a safe way and in accordance with national and local regulations.
7. The Stop Solution contains 0.5 mol/L sulfuric acid and can cause irritation or burns to the skin and eyes. If contact occurs, rinse immediately with plenty of water and seek medical advice.
8. Do not interchange components from kits with different batch numbers. The components have been standardized as a unit for a given batch.
9. Hemolyzed, hyperlipemic, heat-treated or

## STABILITY AND STORAGE

1. Store the kit with all reagents at 2-8°C. Do not freeze.
2. Use all reagents before the expiry date on the kit box label.
3. Diluted Wash Solution Concentrate remains stable for 4 weeks at 2-8°C. If not all wells are to be used, dilute only the portion of Wash Solution Concentrate required.
4. Diluted Sample Diluent Concentrate remains stable for 24 hours at 2-8°C. If not all wells are to be used, dilute only the portion of Sample Diluent Concentrate required.
5. For subsequent use, store unused wells in the foil pouch with the desiccant provided and reseal. Always allow foil pouch to equilibrate to room temperature before opening to avoid condensation in/on the coated microwells.

## COLLECTION OF SPECIMENS

**Handle and dispose of all blood-derived or urine specimens as if they were potentially infectious. See Precautions, sections 2, 4 and 5.**

Determination of NGAL in a single specimen requires 10 µL of urine or plasma. For the assessment of acute kidney injury, serum specimens are not recommended, as NGAL is released from neutrophils during blood clotting and may augment values irrespective of kidney injury<sup>23</sup>. Blood specimens should be collected aseptically

into EDTA or heparin tubes by qualified staff using approved venipuncture techniques. Plasma should be prepared by standard techniques for clinical laboratory testing. Urine should be centrifuged. Cap the prepared clinical specimens. If the assay cannot be performed within 24 hours or specimens are to be shipped, freeze the specimens at  $-20^{\circ}\text{C}$  or below. For long-term storage of clinical specimens,  $-70^{\circ}\text{C}$  or below is recommended. Avoid repeated freezing and thawing. Do not use hemolyzed, hyperlipemic, heat-treated or contaminated specimens.

### PREPARATION OF REAGENTS

1. Bring all specimens and reagents to room temperature ( $20\text{-}25^{\circ}\text{C}$ ). Mix specimens thoroughly by gentle inversion and if necessary clear visible particulate matter by low-speed centrifugation.
2. Determine the number of specimens to be tested (in duplicate) plus any internal laboratory control specimens (in duplicate) plus any reagent blank wells. The pre-coated wells can be used as strips of 8 or as individual wells. Single wells are handled by breaking individual wells apart and placing each well in the frame at an appropriate position. Letters and notches on the wells allow the individual wells to be identified. Add 12 wells for the 6 calibrators (in duplicate). Remove the number of microwells required and replace the remainder in the foil pouch with desiccant at  $2\text{-}8^{\circ}\text{C}$ .
3. Wash Solution: Dilute the 25x Wash Solution Concentrate by pouring the total contents of the bottle (30 mL) into a 1-L graduated cylinder and add distilled or deionized water to a final volume of 750 mL. Mix thoroughly and store at  $2\text{-}8^{\circ}\text{C}$  after use. If not all the wells are to be used, dilute only the required volume of Wash Solution Concentrate with 24 volumes of water to produce a 1/25 dilution.
4. Sample Diluent: Dilute the 5x Sample Diluent Concentrate (contains red dye to aid pipetting) by pouring the total contents of the bottle (60 mL) into a 500-mL graduated cylinder and add distilled or deionized water to a final volume of 300 mL. Mix thoroughly and store at  $2\text{-}8^{\circ}\text{C}$  after use. If not all the wells are to be used, dilute only the required volume of Sample Diluent Concentrate with 4 volumes of water to produce a 1/5 dilution.
5. NGAL Rapid Calibrators (contains red dye to aid pipetting): The assigned concentrations are indicated on their labels. Do not dilute further.
6. HRP-conjugated NGAL Antibody (ready to use): Do not dilute further.
7. TMB Substrate (ready to use): Do not dilute further.
8. Stop Solution (ready to use): Do not dilute further.

### PREPARATION OF SAMPLES

The range of the standard curve is 0.2-20 ng/mL and the diagnostically relevant range is 100-500 ng/mL for plasma and 50-500 ng/mL for urine (see Interpretation of Results). Therefore an initial screening dilution of 1/100 for plasma and 1/50 for urine is recommended. The 1/100 dilution can be prepared by diluting 10  $\mu\text{L}$  of plasma in 990  $\mu\text{L}$  of Sample Diluent and the 1/50 dilution can be prepared by diluting 10  $\mu\text{L}$  of urine in 490  $\mu\text{L}$  of Sample Diluent. Dilutions are mixed by inversion or moderate vortexing. Re-assay of out-of-range samples at lower or higher dilution may be necessary. Dilutions lower than 1/25 should not be used.

### ASSAY PROCEDURE

1. Prepare the assay protocol, assigning the appropriate wells for setting up calibrators, diluted patient specimens and any internal laboratory controls in duplicate. If a reference wavelength of 650 or 620 nm is not available on the ELISA reader, a reagent blank well can be assigned. This is set up with 50  $\mu\text{L}$  of Sample Diluent instead of diluted sample and processed like the other wells.
2. Dilute samples according to the expected NGAL concentrations (1/100 for plasma and 1/50 for urine will be suitable for most samples).
3. Pipette a sufficient volume of each calibrator,


each diluted sample and any internal laboratory controls into the appropriate wells of the polypropylene U-microwell plate to permit subsequent transfer of 50  $\mu\text{L}$  volumes to corresponding coated microwells.

4. Pipette 50  $\mu\text{L}$  volumes of HRP-conjugated NGAL Antibody into the corresponding positions in the coated microwells. Then with a multichannel pipette rapidly transfer 50  $\mu\text{L}$  volumes of the calibrator solutions, diluted samples and internal controls from the U-wells into the corresponding coated wells already containing the detection antibody. This method of sample addition is recommended to reduce the difference in incubation time between the first and last samples added to the coated microwells.
5. Cover the wells and incubate for **30 minutes** at room temperature on a shaking platform set at 200/minute.
6. Aspirate the contents of the microwells and wash the microwells three times with at least 300  $\mu\text{L}$  of the previously diluted Wash Solution. If washing is performed manually, empty the microwells by inversion and gentle shaking into a suitable container, followed by blotting in the inverted position on a paper towel. A dwell time of 1 minute before emptying is recommended for at least the last wash of the cycle. The vigor with which diluted Wash Solution is filled into or emptied from the wells influences final color development. Manual pipetting, which may be very gentle and lead to high color development, is only recommended in the absence of alternatives such as filling the wells by immersion, using a multichannel manual washing dispenser, or using an automatic washing apparatus.
7. Dispense 100  $\mu\text{L}$  of TMB Substrate (ready to use) into each microwell. The use of a multichannel micropipette is recommended to reduce pipetting time. Cover the wells and incubate for **exactly 15 minutes** at room temperature in the dark. Start the clock when filling the first well.
8. Add 100  $\mu\text{L}$  of Stop Solution (ready to use)

to each well, maintaining the same pipetting sequence and rate as in Step 7. Mix by gentle shaking for 20 seconds, avoiding splashing. Read the wells within 30 minutes.

9. Read the absorbances of the wells at 450 nm in an appropriate microplate reader (reference wavelength 650 or 620 nm). If no reference wavelength is available, the value of the reagent blank well is subtracted from each of the other values before other calculations are performed.

#### SCHEMATIC OVERVIEW

## CALCULATION OF RESULTS

A calibration curve is constructed by plotting the mean of duplicate absorbance values for each NGAL Rapid Calibrator on the y-axis against the corresponding NGAL concentrations in ng/mL on the x-axis. The calibration curve must meet the validation requirements. The NGAL concentration of each diluted sample is then found by locating the point on the curve corresponding to the mean of duplicate absorbance values for the diluted sample and reading its corresponding concentration in ng/mL from the x-axis. The concentration of NGAL in the undiluted specimen is calculated by multiplying this result by the sample dilution factor.

This procedure can be performed manually using graph paper with linear x and y scales. A smooth curve can be drawn through the points or adjacent points can be joined by straight lines. The latter procedure may slightly overestimate concentration values between points when the curve is slightly convex to the left, which is the typical finding. Although the curve may approximate to a straight line, it is both practically and theoretically incorrect to calculate and draw the straight line of best fit and to read the results from this.

The procedure can also be performed by an ELISA reader software program incorporating curve fitting procedures. The procedure of choice is to use linear x and y axes with 4-parameter logistic curve fitting. Diluted samples that give a mean absorbance above that for the NGAL Rapid Calibrator 6 or below that for the NGAL Rapid Calibrator 2 are out of the range of the assay and their concentrations should be noted as  $>20$  ng/mL and  $<0.2$  ng/mL, respectively. The corresponding concentrations in the undiluted samples are calculated  $>(20 \times \text{dilution factor})$  ng/mL and  $<(0.2 \times \text{dilution factor})$  ng/mL, respectively. If necessary, these samples can be re-assayed at higher and lower dilutions for high- and low-reading samples, respectively. The new dilution factors should be those estimated to give absorbance values that fall well within the range of the calibration curve, but dilutions lower than 1/25 should not be used.

## VALIDATION OF CALIBRATION CURVE

The mean absorbance for the NGAL Rapid Calibrator 6 should be  $>1.5$ . The mean absorbance for any NGAL Rapid Calibrator should be higher than that for the previous NGAL Rapid Calibrator, e.g. absorbance(NGAL Rapid Calibrator 6)  $>$  absorbance(NGAL Rapid Calibrator 5). The curve should be slightly convex to the left when the results are plotted on linear axes.

## CALIBRATION TROUBLESHOOTING

**High calibration curve:** In some circumstances (e.g. high ambient temperature, gentle washing technique) the absorbance values of the calibration curve may be generally elevated, so that the curve flattens off at the higher concentrations giving values near the upper reading range of the microplate reader. This can be compensated for by re-reading all the wells at 405 nm, giving rise to lower absorbance values within the more linear response range of the reader. Results are then calculated from the readings made at the 405 nm wavelength. This avoids having to repeat the assay, but because wells should be read within 30 minutes of adding Stop Solution, it is recommended that the initial readings obtained at 450 nm should be examined within this time period to see whether it is necessary to re-read the plate at 405 nm. If it is known that local laboratory conditions regularly give rise to high calibration curves, the problem can also be overcome by shortening the color development time to 10 or even 5 minutes.

## Out-of-line points for individual calibrators:

One or more individual calibrators may give anomalous absorbance readings. One or both of the duplicate values may be out of line, and the mean of the duplicates may be out of line. This error is significant if it impairs satisfactory curve fitting by the 4-parameter logistic method, which, as a result of the anomalous value, is shifted away from other calibrator points that are in fact correct. The calibrator points and fitted curve should always be examined for correct fit before any calculations of

concentration from it are accepted. A poorly fitting curve will also be revealed by a high sum of residual squares. If only one calibrator is affected, which is not the highest calibrator, two courses of action are possible:

i) An erroneous singlet or duplicate result should be eliminated from the curve, and the remaining results refitted by the 4-parameter logistic procedure. If a satisfactory fit is obtained, provisional concentration results can be calculated from it.

ii) If no satisfactory fit can be obtained in this way, but the curve is otherwise consistent, provisional results can be obtained from straight lines or simple cubic spline fitting between the means of duplicates, omitting the erroneous point.

If two or more calibrators are affected, the assay should be repeated.

A deviant result for an individual calibrator can be due to operator error or to calibrator deterioration. If both duplicate values are consistently out of line in successive assays, the calibrator is faulty and should be omitted.

#### TRACEABILITY OF CALIBRATOR VALUE

No internationally approved reference material for NGAL is currently available. The NGAL concentration of calibrator material has been assigned by turbidimetry using a precise transfer protocol ensuring traceability to the BioPorto Diagnostics master calibrator. The reference material was value-assigned by measurement of light absorbance at 280 nm using a theoretically calculated extinction coefficient based on the amino-acid composition.

#### INTERPRETATION OF RESULTS

The finding of a raised urinary or plasma level of NGAL cannot be independently diagnostic of any single pathology. As stated in the Introduction, a variety of independent pathologies are associated with raised levels of urinary or plasma NGAL. Physicians must interpret the significance of any raised NGAL level in the light of each patient's clinical features.

Urinary and/or plasma concentrations of NGAL may be substantially elevated in conditions that

have not affected the kidney, including bacterial infections, other inflammatory disorders and certain carcinomas. BioPorto Diagnostics has determined that the NGAL concentration in an isolated sample of urine or plasma should exceed 250 ng/mL in order to indicate the presence of renal disorder, including acute kidney injury, without incurring the risk of an unacceptably high proportion of false positive diagnoses of renal disorder.

#### QUALITY CONTROL

Laboratories intending to perform repeated assays should establish their own high-reading and low-reading control samples, stored in small (e.g. 50 µL) aliquots at -70°C or below. An aliquot of each should be thawed and tested in each assay and a record kept of successive results. This serves as a control of test performance, test integrity and operator reliability. The results should be examined for drift (tendency for successive results to rise or fall) or significant deviation from the mean of previous results. Values not deviating by more than 20% from the mean of previous values can be taken to indicate acceptability of the assay. Aliquots of controls should not be refrozen for repeated assay once thawed, and if a further assay is performed, fresh control aliquots and fresh dilutions of patient specimens should be used.

#### LIMITATIONS

The finding of a raised urinary or plasma level of NGAL cannot be independently diagnostic of any single pathology. As stated in the Introduction, a variety of independent pathologies are associated with raised levels of urinary or plasma NGAL. Physicians must interpret the significance of any raised NGAL level in the light of each patient's clinical features.

#### EXPECTED RESULTS

The mean NGAL concentration in samples from healthy donors was 63 ng/mL (range 37–106 ng/mL, n = 80) in EDTA plasma and 5.3 ng/mL (range 0.7–9.8 ng/mL, n = 7) in urine. In unselected patients admitted to intensive care, the NGAL concentrations in urine ranged from 9 ng/mL to 40,000 ng/mL

(40 µg/mL) in urine (n = 60) and from 25 ng/mL to 3490 ng/mL in EDTA plasma (n = 60).

### PERFORMANCE CHARACTERISTICS

**Limit of detection:** The lowest concentration of NGAL giving an absorbance reading greater than 2 SD above the mean zero (NGAL Rapid Calibrator 1) reading (n = 8) was determined to be 0.008 ng/mL, this being lower than the value of NGAL Rapid Calibrator 2.

**Precision:** Intraassay variation was determined by measurement of NGAL in two urine samples and two plasma samples with 8 replicates. The following results were obtained (CV = coefficient of variation):

Samples	CV
Urine A	3.4%
Urine B	4.3%
Plasma A	2.9%
Plasma B	1.9%

Interassay variation was determined by measurement of NGAL in 2 diluted urine samples and 2 diluted EDTA plasma samples with 2 replicates in 4 separate assays. The following results were obtained:

Samples	CV
Urine A	4.7%
Urine B	22.7%
Plasma A	11.4%
Plasma B	12.4%

**Analytical recovery:** Urine and plasma samples were spiked with recombinant human NGAL and analyzed in the assay. Recovery was calculated from (Measured/Expected) expressed as a percentage.

Sample	Measured (range)	Expected (range)	Recovery (range)
Urine	4.3-12.0 ng/mL	4.5-12.4 ng/mL	93%-98%
Plasma	4.2-10.5 ng/mL	4.6-12.6 ng/mL	83%-90%

**Linearity:** NGAL was measured in serial dilutions (n = 8) of 2 urine samples and 2 plasma samples. The CV of the mean of the measured values corrected for the dilution was 9.9% and 12.5% for the two urine samples and 12.2% and 7.4% for the two plasma samples, demonstrating satisfactory linearity.

**Sample material:** Analysis of samples of urine or EDTA plasma (or serum) showed no significant differences in analytical recovery, linearity or precision. However, serum samples are not recommended for the assessment of acute kidney injury as NGAL is released from neutrophils during blood clotting and may augment values irrespective of kidney injury<sup>23</sup>.

**Specificity:** The two monoclonal antibodies against human NGAL used in the assay have been shown to bind to different preparations of recombinant human NGAL and to give a single band at 25 kDa on Western blot analysis of a reduced post-nuclear supernatant from human neutrophils<sup>25</sup>.

### LIABILITY

This ELISA Kit is only intended for the in vitro determination of NGAL in human urine or plasma. The ELISA Kit is only intended for use by qualified personnel carrying out research or diagnostic activities. If the recipient of this test passes it on in any way to a third party, this instruction must be enclosed, and said recipient shall at recipient's own risk secure in favor of BioPorto Diagnostics A/S all limitations of liability herein.

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1. Allen RA, Erickson RW, Jesaitis AJ (1989) Identification of a human neutrophil protein of Mr 24000 that binds *N*-formyl peptides: co-sedimentation with specific granules. *Biochim Biophys Acta* 991:123-133.
2. Kjeldsen L, Johnsen AH, Sengelov H, Borregaard N (1993) Isolation and primary structure of NGAL, a novel protein associated with human neutrophil gelatinase. *J Biol Chem* 268:10425-10432.
3. Xu SY, Carlson M, Engstrom A, Garcia R, Peterson CG, Venge P (1994) Purification and characterization of a human neutrophil lipocalin (HNL) from the secondary granules of human neutrophils. *Scand J Clin Lab Invest* 54:365-376.
4. Flower DR, North AC, Attwood TK (1991) Mouse oncogene protein 24p3 is a member of the lipocalin protein family. *Biochem Biophys Res Commun* 180:69-74.
5. Liu Q, Ryon J, Nilsen-Hamilton M (1997) Uterocalin: a mouse acute phase protein expressed in the uterus around birth. *Mol Reprod Dev* 46:507-514.
6. Stoesz SP, Gould MN (1995) Overexpression of neu-related lipocalin (NRL) in neu-initiated but not *ras* or chemically initiated rat mammary carcinomas. *Oncogene* 11:2233-2241.
7. Triebel S, Blaser J, Reinke H, Tschesche H (1992) A 25 kDa alpha 2-microglobulin-related protein is a component of the 125 kDa form of human gelatinase. *FEBS Lett* 314:386-388.
8. Yan L, Borregaard N, Kjeldsen L, Moses MA (2001) The high molecular weight urinary matrix metalloproteinase (MMP) activity is a complex of gelatinase B/MMP-9 and neutrophil gelatinase-associated lipocalin (NGAL). Modulation of MMP-9 activity by NGAL. *J Biol Chem* 276:37258-37265.
9. Cowland JB, Borregaard N (1997) Molecular characterization and pattern of tissue expression of the gene for neutrophil gelatinase-associated lipocalin from humans. *Genomics* 45:17-23.
10. Friedl A, Stoesz SP, Buckley P, Gould MN (1999) Neutrophil gelatinase-associated lipocalin in normal and neoplastic human tissues. Cell type-specific pattern of expression. *Histochem J* 31:433-441.
11. Nielsen BS, Borregaard N, Bundgaard JR, Timshel S, Sehested M, Kjeldsen L (1996) Induction of NGAL synthesis in epithelial cells of human colorectal neoplasia and inflammatory bowel diseases. *Gut* 38:414-420.
12. Stoesz SP, Friedl A, Haag JD, Lindstrom MJ, Clark GM, Gould MN (1998) Heterogeneous expression of the lipocalin NGAL in primary breast cancers. *Int J Cancer* 79:565-572.
13. Monier F, Surla A, Guillot M, Morel F (2000) Gelatinase isoforms in urine from bladder cancer patients. *Clin Chim Acta* 299:11-23.
14. Bangert K, Ghiglionone M, Heslet L, Uttenhalt LO (2005) Urinary NGAL is dramatically increased in acute renal failure. *Intensive Care Med* 31(Suppl 1):S65.
15. Goetz DH, Holmes MA, Borregaard N, Bluhm ME, Raymond KN, Strong RK (2002) The neutrophil lipocalin NGAL is a bacteriostatic agent that interferes with siderophore-mediated iron acquisition. *Mol Cell* 10:1033-1043.
16. Tong Z, Wu X, Ovcharenko D, Zhu J, Chen CS, Kehrher JP (2005) Neutrophil gelatinase associated lipocalin as a survival factor. *Biochem J* 391:441-448.
17. Hraba-Renevey S, Turler H, Kress M, Salomon C, Weil R (1989) SV40-induced expression of mouse gene 24p3 involves a post-transcriptional mechanism. *Oncogene* 4:601-608.
18. Ohlsson S, Wieslander J, Segelmark M (2003) Increased circulating levels of proteinase 3 in patients with anti-neutrophilic cytoplasmic auto-antibodies-associated systemic vasculitis

- in remission. *Clin Exp Immunol* 131:528-535.
19. Xin C, Yulong X, Yu C, Changchun C, Feng Z, Xinwei M (2008) Urine neutrophil gelatinase-associated lipocalin and interleukin-18 predict acute kidney injury after cardiac surgery. *Renal fail* 30:904-913.
  20. Tuladhar SM, Püntmann VO, Soni M, Punjabi PP, Bogle RG (2009) Rapid detection of acute kidney injury by plasma and urinary neutrophil gelatinase-associated lipocalin after cardiopulmonary bypass. *J Cardiovasc Pharmacol* 53:261-266.
  21. Bangert K, Heslet L, Ghigliione M, Uttenthal LO (2006) NGAL is significantly increased in urine and plasma in acute renal failure. *Intensive Care Med* 32(Suppl 1):S10.
  22. Bolignano D, Donato V, Coppolino G, Campo S, Buemi A, Lacquaniti A, Buemi M (2008) Neutrophil gelatinase-associated lipocalin (NGAL) as a marker of kidney damage. *Am J Kidney Dis* 52:595-605.
  23. Xu SY, Pauksen K, Venge P (1995) Serum measurements of human neutrophil lipocalin (HNL) discriminate between acute bacterial and viral infections. *Scand J Clin Lab Invest* 55:125-131.
  24. Fernandez CA, Yan L, Louis G, Yang J, Kutok JL, Moses MA (2005) The matrix metalloproteinase-9/neutrophil gelatinase-associated lipocalin complex plays a role in breast tumor growth and is present in the urine of breast cancer patients. *Clin Cancer Res* 11:5390-5395.
  25. Kjeldsen L, Koch C, Arnljots K, Borregaard N (1996) Characterization of two ELISAs for NGAL, a newly described lipocalin in human neutrophils. *J Immunol Methods* 196:155-164.

Catalogue number

In vitro diagnostic medical device

Batch code

Consult instructions for use

European Conformity

Use by

Manufacturer

Keep away from sunlight

Temperature limitation

Do not reuse

Caution, consult accompanying documents

Biological risk

Do not use if package is damaged

Concentrated Wash Solution.  
Dilute before use.

Concentrated Sample  
Diluent. Dilute before use.

**BIOPORTO**<sup>®</sup>  
Diagnostics

BioPorto Diagnostics A/S  
Tuborg Havnevej 15, st.  
DK-2900 Hellerup  
Denmark

Phone (+45) 4529 0000  
Fax (+45) 4529 0001  
E-mail [info@bioporto.com](mailto:info@bioporto.com)  
Web [www.bioporto.com](http://www.bioporto.com)  
[www.ngal.com](http://www.ngal.com)