

APPLICATION NOTE FOR

ROCHE MODULAR® P¹

The NGAL Test™ Reagent Kit

REF/Cat.No.	ST001RA		ST002RA	ST003RA
Product name	The NGAL Test™ Reagent Kit		The NGAL Test™ Calibrator Kit	The NGAL Test™ Control Kit
	R1	R2	50, 150, 600, 1500, 3000 ng/mL	Low and High
	1 x 35 mL	1 x 7 mL	5 x 1 mL	3 x 1 mL x 2 levels

Number of determinations: 1 mL of immunoparticle suspension **R2** provides 20 cuvette readings with the provided settings in this application. The dead volume of the analyzer and reagent container should be added when calculating the required amount of reagent.

PERFORMANCE DATA

The performance data shown were obtained by the manufacturer for this particular analyzer model. For additional performance data and product application, please read the instructions for use accompanying the product carefully. Each individual laboratory should validate the use of The NGAL Test™ on its system.



CALIBRATION STABILITY

It is recommended to recalibrate every 4 weeks, when reagent lots change or quality control results fall outside the range as established by the individual laboratory.

TROUBLE SHOOTING

If performance is unacceptable, try to recalibrate. Check reagents and procedure. If the problem persists, please contact instrument supplier or reagent supplier.

¹ Modular® P is a registered trademark of Roche Diagnostics GmbH, Mannheim, Germany.

APPLICATION PARAMETERS

Assay/Time/Point	2-Point-End	10	18	34	0	0			
Wave (2nd/Primary)	800	570							
Sample Volume			Reagent Vol.						
Normal	3.0	0.0	0	R1	150	0	xxx	0	Timing
Decrease	15.0	3.0	105	R2	0	0	xxx	0	T2
Increase	6.0	0.0	0	R3	50	0	xxx	0	T3
				R4	0	0	xxx	0	
Diluent									
<input type="radio"/> Water <input checked="" type="radio"/> Diluent			314	0					
Abs. Limit	32000	Increase							
Prozone Limit	0	0	0	0	0	0	Lower		
Cell Detergent	Detergent 1								
Twin Test	Cancel								

CALIBRATION		Auto Calibration																						
<table border="1"> <tr> <td>Calibration Type</td> <td>SPLINE</td> </tr> <tr> <td>Point</td> <td>6</td> </tr> <tr> <td>Span</td> <td>6</td> </tr> <tr> <td>Weight</td> <td>0</td> </tr> <tr> <td>Update Type</td> <td>none</td> </tr> <tr> <td>Isozyme Q Channel</td> <td>Cancel</td> </tr> </table>		Calibration Type	SPLINE	Point	6	Span	6	Weight	0	Update Type	none	Isozyme Q Channel	Cancel	<table border="1"> <tr> <td>Change Over</td> <td>Cancel</td> </tr> <tr> <td>Modul</td> <td>Cancel</td> </tr> <tr> <td>Lot</td> <td>Cancel</td> </tr> <tr> <td>Bottel</td> <td>Cancel</td> </tr> </table>		Change Over	Cancel	Modul	Cancel	Lot	Cancel	Bottel	Cancel	
Calibration Type	SPLINE																							
Point	6																							
Span	6																							
Weight	0																							
Update Type	none																							
Isozyme Q Channel	Cancel																							
Change Over	Cancel																							
Modul	Cancel																							
Lot	Cancel																							
Bottel	Cancel																							
<table border="1"> <tr> <td>SD Limit</td> <td>999</td> </tr> <tr> <td>Duplicate Limit</td> <td>99</td> <td>32000</td> </tr> <tr> <td>Sensitivity Limit</td> <td>-99999</td> <td>99999</td> </tr> <tr> <td>S1-Abs Limit</td> <td>-32000</td> <td>32000</td> </tr> </table>		SD Limit	999	Duplicate Limit	99	32000	Sensitivity Limit	-99999	99999	S1-Abs Limit	-32000	32000	<table border="1"> <tr> <td colspan="2">Timeout</td> </tr> <tr> <td><input type="checkbox"/> Blank</td> <td>1</td> <td>2 Point</td> <td>1</td> </tr> <tr> <td><input type="checkbox"/> Span</td> <td>1</td> <td>Full</td> <td>1</td> </tr> </table>		Timeout		<input type="checkbox"/> Blank	1	2 Point	1	<input type="checkbox"/> Span	1	Full	1
SD Limit	999																							
Duplicate Limit	99	32000																						
Sensitivity Limit	-99999	99999																						
S1-Abs Limit	-32000	32000																						
Timeout																								
<input type="checkbox"/> Blank	1	2 Point	1																					
<input type="checkbox"/> Span	1	Full	1																					
<table border="1"> <tr> <td colspan="2">V Auto Masking</td> </tr> <tr> <td>QC Violation</td> <td></td> </tr> <tr> <td>Method</td> <td>Blank</td> <td>Rule</td> <td>1s</td> </tr> <tr> <td>Control</td> <td>None</td> <td></td> <td></td> </tr> <tr> <td></td> <td>None</td> <td></td> <td></td> </tr> <tr> <td></td> <td>None</td> <td></td> <td></td> </tr> </table>		V Auto Masking		QC Violation		Method	Blank	Rule	1s	Control	None				None				None					
V Auto Masking																								
QC Violation																								
Method	Blank	Rule	1s																					
Control	None																							
	None																							
	None																							

RANGE

Application Code xxx
 Unit ng/mL

Report Name NGAL
 Data Mode Active
 V Automatic rerun _____
 Technical Limit 25 3000
 Repeat Limit -99999 999999
 Control Interval Time 1

Qualitative Control Interval Time

(1)	0	--
(2)	0	-
(3)	0	+/-
(4)	0	++
(5)	0	+++
(6)		++++

Expected Ranges

Male

<u>99</u>	Year	<u>-99999</u>	<u>99999</u>
<u>100</u>	Year	<u>-99999</u>	<u>99999</u>
		<u>-99999</u>	<u>99999</u>

Female

<u>99</u>	Year	<u>-99999</u>	<u>99999</u>
<u>100</u>	Year	<u>-99999</u>	<u>99999</u>
		<u>-99999</u>	<u>99999</u>

Default

Sex Male Female

Range Range 1 Range 2 Range 3

OTHERS

Standards

	(1)	(2)	(3)	(4)	(5)	(6)
Calibrator Code	<u>501</u>	<u>xxx</u>	<u>xxx</u>	<u>xxx</u>	<u>xxx</u>	<u>xxx</u>
Concentration	<u>0</u>	<u>50</u>	<u>150</u>	<u>600</u>	<u>1500</u>	<u>3000</u>
Rack No. – Pos.	<u>XXXXX</u>	<u>XXXXXX</u>	<u>XXXXXX</u>	<u>XXXXXX</u>	<u>XXXXXX</u>	<u>XXXXXX</u>
Sample Volume	<u>3.0</u>	<u>3.0</u>	<u>3.0</u>	<u>3.0</u>	<u>3.0</u>	<u>3.0</u>
Diluent S. Volume	<u>0.0</u>	<u>0.0</u>	<u>0.0</u>	<u>0.0</u>	<u>0.0</u>	<u>0.0</u>
Diluent Volume	<u>0</u>	<u>0</u>	<u>0</u>	<u>0</u>	<u>0</u>	<u>0</u>

xxx: To be defined by operator.