For your clinical chemistry analyzer

PERFORMANCE DATA AND APPLICATION NOTE FOR SIEMENS ADVIA® 2400¹

The NGAL Test™ Reagent Kit

REF/Cat. No.	ST001RA	ST002RA	ST003RA
Product name	The NGAL Test™ Reagent Kit	The NGAL Test™ Calibrator Kit	The NGAL Test™ Control Kit
	R1	R2	
	150, 600, 1500, 3000, 5000 ng/mL	150, 600, 1500, 3000, 5000 ng/mL	Low and High
	1 x 35 mL	1 x 7 mL	5 x 1 mL
			3 x 1 mL x 2 levels

Number of determinations: 1 mL of immunoparticle suspension **R2** provides 50 cuvette readings with the provided settings in this application. The dead volume of the analyzer and reagent container should be added when calculating the required amount of reagent.

Regulatory status:
For Research Use Only. Not for use in diagnostic procedures.

CALIBRATION

The NGAL Test™ Calibrator Kit (**REF** ST002RA) should be used the following way:

Calibration point	Kit calibrator	Level
1	Use saline as blank*	0 ng/ml
2	Calibrator 1 (150 ng/mL)**	50 ng/mL
3	Calibrator 1 (150 ng/mL)	150 ng/mL
4	Calibrator 2 (600 ng/mL)	600 ng/mL
5	Calibrator 3 (1500 ng/mL)	1500 ng/mL
6	Calibrator 4 (3000 ng/mL)	3000 ng/mL

* Saline is not included in the kit

** Special dilution, consult settings below

NB: ST002RA Calibrator 5 (5000 ng/mL) is not used for calibration on Siemens ADVIA® 2400 system.

CALIBRATION STABILITY

It is recommended to recalibrate every 4 weeks, when reagent lots change or quality control results fall outside the range as established by the individual laboratory.

TROUBLE SHOOTING

If performance is unacceptable, try to recalibrate. Check reagents and procedure. If the problem persists, please contact instrument supplier or reagent supplier.

1. ADVIA® is a registered trademark of Siemens Healthcare Diagnostics Inc., New York, USA

The NGAL Test™

For your clinical chemistry analyzer

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* Use saline as blank solution

** Use ST002RA Calibrator 1 (150 ng/mL) to make the 1st calibration point 50 ng/mL

*** ST002RA Calibrator 5 (5000 ng/mL) is not used for calibration on Siemens ADVIA® 2400 system